

LE 820 Circular Pool

Provides the ideal environment to carry out the Swimming Test described by Morris (Morris .G.M., Learning and Motivation, 12:239, 1981)



The island can be located anywhere in the pool. With the pool are supplied two platforms, 80 and 110 mm diameter.

The pool are manufactured of Polypropylene. It's located in a support with 4 wheels for better displacement.

MODELS

<u>Ref</u>	<u>Diam. (cm)</u>	<u>High (cm)</u>
LE 820/90	90	50
LE 820/120	120	60
LE 820/140	140	60
LE 820/170	170	60
LE 820/200	200	60



In the Control Box are located:

- The Heater
- The Water Circulation Pump
- The Level Controller
- Electro valve for Pool filling

The Level Controller acts directly on the Electro Valve, turning it off when the liquid arrives to the corresponding high.

TECHNICAL CHARACTERISTICS

Heater Power:	3.000 W
Temperature Regulation:	22 - 32°C (depending on environment conditions)
Heating Speed:	3°C / h (in the biggest model)
Power Supply:	220 V, 50 Hz
input / Output Intake:	19 mm